

M-NEB® VET INHALATION SET FOR SMALL ANIMALS

The M-neb® vet inhalation set is used to treat respiratory diseases in small animals such as small dogs, cats, rabbits, hedgehogs, falcons and reptiles. With the help of innovative mesh technology, it enables **targeted and effective** administration of a wide range of inhalation solutions, which reach the **deepest parts of your pet's airways** through inhalation.

In addition to the centrepiece, the control and mesh nebulizer unit, the set contains an inhalation box and a small inhalation mask. The most suitable type of application for your pet can be selected according to your **individual needs**.

The M-neb® vet inhalation set is also **ideal for anxious animals** thanks to its virtually silent and very pleasant operation.

NEBU-TEC: Innovative inhalation systems for humans and animals

NEBU-TEC is a leading, internationally active medical technology company based in Elsenfeld (Germany). With around 40 dedicated employees, we develop, produce and distribute modern inhalation systems based on advanced mesh technology. The inhalation systems are used in human and veterinary medicine. 95% of the production takes place in Germany – NEBU-TEC therefore stands for the highest quality and reliability.

NEBU-TEC is considered a pioneer in veterinary medicine and offers specialized inhalation systems for the treatment of horses, small animals such as dogs, cats and falcons.

CONTACT

Would you like to support your pet with an effective inhalation device? Contact us by phone or e-mail – **we will be happy to advise you:**



NEBU-TEC International med. Produkte Eike Kern GmbH

Kreuzfeldring 17 | D - 63820 Elsenfeld

+49 (0) 6022 - 610 62 -0 | info@nebu-tec.de

www.nebu-tec.de

You can also visit us on
our social media channels:



NTL_HM-300_Productfolder-EN_LB | 01.02.2025



STRESS-FREE INHALATION THERAPY FOR SMALL DOGS, CATS AND OTHER SMALL ANIMALS



APPLICATION & USE OF INHALATION

Respiratory diseases are among the **most frequently diagnosed** diseases in veterinary medicine. The M-neb® vet inhalation set can be used for a wide variety of small animals, for example ...

... for small dogs

The most common respiratory disease in dogs is **chronic bronchitis**, which usually develops from **acute bronchitis** that has not been (sufficiently) treated. Dogs also suffer from normal colds, but **kennel cough, pneumonia and allergic asthma** are also common.

... for cats

Respiratory diseases are relatively common in cats and can occur at any age. The causes vary depending on the disease and range from an infection with **bacteria or viruses, allergies** and exposure to pollutants to **poor housing conditions**. Cat owners are familiar with the so-called cat cold complex, but also diseases of the lower respiratory tract, such as **bronchitis, pneumonia or feline asthma**.

... for other small animals

Other small animals, such as **rabbits or hedgehogs, as well as falcons**, also suffer from respiratory diseases that may require inhalation treatment. **Respiratory symptoms** are also a common reason for **reptiles** to be seen in veterinary practices.



FUNCTIONALITY

The form of application selected according to the animal species and **individual requirements** is supplemented by the use of the M-neb® vet mesh nebulizer.

M-neb® vet: smart nebulization at the touch of a button

The leak-proof **mesh nebulizer** unit filled with the inhalation solution is connected to either the **inhalation box** or the **inhalation mask** and the **M-neb® vet control unit** using the supplied connection cable. The control unit is operated by the pet owner and has two modes, **intermittent and continuous**, which can be selected according to the inhalation solution used. Once the solution has been completely nebulized, **the device switches off automatically**.

The M-neb® vet mesh nebulizer unit produces **optimal particle sizes for treating the upper and lower respiratory tract of small animals**. The entire inhalation system works **almost silently**, so that even anxious pets can get used to inhalation therapy **quickly and easily**.



THE WAY YOUR PET LOVES IT MOST ...

The M-neb® vet inhalation set can be used depending on the behaviour and character of the small animal to be treated and thus according to the individual needs with the application form of your choice:

Inhalation via a small inhalation mask

If you want to give your pet a sense of presence, the **small inhalation mask** is ideal for you. You can hold your pet for example in your **arms or on your lap**. You control the inhalation session manually and are therefore able to respond individually to your pet's **breathing rhythm**. The application directly in front of the mouth area is **gentle on the eyes and skin**, so this method has proven effective for **administering medication**. Direct inhalation also **prevents loss of active ingredients**.

Also for anxious small animals

Silent and very pleasant operation leads to rapid acceptance

From home without much effort

Easy to use as well as quick and hygienic cleaning for long-term and effective inhalation therapy



Inhalation in the practical inhalation box

The inhalation box is specially designed for pets that are **hesitant in new situations** and prefer not to be held. Your pet remains **relaxed in the inhalation box** throughout the entire therapy session, breathing in the aerosol that collects inside. Because this method not only treats the respiratory tract but also the largest respiratory organ – the skin – through **room nebulization**, this application is particularly suitable for the **nebulization of saline solutions**.

Precise and accurate

Inhalation via two different modes focussing on your pet's breathing rhythm to achieve the best results

Many preparations can be nebulized

All common inhalation solutions can be used without heating and the associated destruction

100 % Made in Germany

We select all the materials we use with great care and only use proven and mostly biocompatible materials from human medicine, which we source from German suppliers for our quality products

TECHNICAL DATA

The M-neb® vet inhalation set for small dogs, cats and other small animals with **modern and sophisticated technology** that will **inspire you**:



L550 x B360 x H330 mm, 2.350 g total net weight (in operational condition)



MMD: > 4 µm | FPF: < 50 %
Output: 0,25 - 0,75 ml per minute (in continuous mode)



Intermittent & continuous mode
Automatic switch-off function
Integrated charge level indicator

