

SADOMA®-II INHALATION SET FOR DOGS

The SaDoMa®-II inhalation set for dogs is a **battery-powered inhalation device** that uses **innovative mesh technology** to convert liquids (NaCl, brine, medication, etc.) into aerosols that reach the deepest parts of your dog's airways.

It consists of the **SaDoMa®-II inhalation mask** specially developed for dogs with a distinct profile and is supplemented by the **M-neb® vet mesh nebulizer**.

The **transparent inhalation mask** is attached to the dog's head with a fastening strap. The aerosol is produced by the mesh nebulizer unit located on the mask directly in front of your dog's snout, reaches its destination directly **without detours or condensation losses**.

With **two different modes**, you can also adapt the nebulization to the **dog's breathing profile** so that your faithful canine companions feels better quickly.

NEBU-TEC: Innovative inhalation systems for humans and animals

NEBU-TEC is a leading, internationally active medical technology company based in Elsenfeld (Germany). With around 40 dedicated employees, we develop, produce and distribute modern inhalation systems based on advanced mesh technology. The inhalation systems are used in human and veterinary medicine. 95% of the production takes place in Germany – **NEBU-TEC** therefore stands for the highest quality and reliability.

NEBU-TEC is considered a pioneer in veterinary medicine and offers specialized inhalation systems for the treatment of horses, small animals such as dogs, cats and falcons.

CONTACT

You would like to support your dog with an effective inhalation device support? Contact us by phone or e-mail – **we will be happy to advise you:**



NEBU-TEC International med. Produkte Eike Kern GmbH

Kreuzfeldring 17 | D - 63820 Elsenfeld

+49 (0) 6022 - 610 62 -0 | info@nebu-tec.de

www.nebu-tec.de

You can also visit us on
our social media channels:



NTL_DM-200_Productfolder-EN_B | 01.02.2025



USE & APPLICATION OF THE INHALATION

Respiratory diseases in dogs are among the **most frequently** diagnosed illnesses. The causes vary and range from an infection with **bacteria** or **viruses**, **allergies** and **exposure to pollutants** to **poor housing conditions**. However, the symptoms are usually very similar. While coughing is usually the main symptom, **increased mucus production**, **breathing problems**, **loss of performance** or **choking** and **wheezing** can also indicate a **respiratory disease**.

The **most common respiratory disease** in dogs is **chronic bronchitis**, which usually develops from acute **bronchitis** that has not been (sufficiently) treated. However, dogs also suffer from normal colds, kennel cough, pneumonia or allergic asthma.

Medication with an anti-inflammatory, bronchodilator or expectorant effect is often used to treat **respiratory diseases** in dogs. As systemic administration in the form of tablets or injections can have serious side effects, especially with cortisone, **inhaled administration of medication** is **much gentler** on the dog.

COMFORTABLE INHALATION FOR CANINE COMPANIONS



FUNCTIONALITY

The SaDoMa®-II inhalation set for dogs consists of the SaDoMa®-II inhalation mask specially developed for dogs with a distinct profile and is supplemented by the M-neb® vet mesh nebulizer.

M-neb® vet: smart nebulization at the touch of a button

The leak-proof mesh nebulizer unit filled with the inhalation solution is connected to both the SaDoMa®-II inhalation mask and the M-neb® vet control unit using the supplied connection cable. The control unit is operated by the dog owner and has two modes, **intermittent and continuous**, which can be selected according to the inhalation solution used (medication, brine, homeopathic remedies, etc.). When the solution is completely nebulized, the **device switches off automatically**.

The M-neb® vet mesh nebulizer unit produces **optimal particle sizes for treating the upper and lower respiratory tract of dogs**. The entire inhalation system works almost silently, so that even anxious dogs can **quickly and easily** become accustomed to inhalation therapy.



Whether for a simple cough as part of a cold or chronic bronchitis - the tried and tested inhalation set developed with veterinary surgeons promises, everything necessary for quick and effective inhalation therapy for dogs:

From home and without much effort

Easy to use as well as quick and hygienic cleaning for long-term and effective inhalation therapy

The latest mesh technology

Optimum particle sizes and excellent output values, adapted to your dog's respiratory tract



Individual Customisation option

The inhalation mask can be easily adapted to different head sizes and dog breeds using different sized sealing rings

Also suitable for anxious dogs

Silent and very pleasant way of working leads to quick acceptance

Ergonomic & transparent

Lightweight inhalation mask with direct aerosol injection directly in front of your dog's snout

100 % Made in Germany

We select all the materials we use with great care and only use proven and mostly biocompatible materials from human medicine, which we source from German suppliers for our quality products

For effortless breathing

Inhalation and exhalation valves allow your dog to breathe easily and prevent losses

TECHNICAL DATA

The SaDoMa®-II inhalation set for dogs is equipped with modern and sophisticated technology that will inspire you:



165 x 135 x 155 mm, 475 g total weight (in operational condition)



MMD: >4 µm | FPF: < 50 %
Output: 0,25 - 0,75 ml per minute (in continuous mode)



Intermittent & continuous mode
Automatic switch-off function
Integrated charge level indicator

